

iCOCOON



Relaxation Chair

HCP-ICO10000A

HCP-ICO10000D

HCP-ICO10000E

Operating Manual

Attached Document

Thank you for purchasing the iCOCOON. This is a high-end, furniture-style relaxation chair that combines design and medication with the goal of promoting mind and body health.

Like a butterfly emerging from its cocoon, sitting in this chair makes you feel reborn into a new self. It's a comfort like no other in the world.

For all customers from 10 to 100 years old, we hope your family will enjoy this chair for many years to come as a necessity that makes you feel happy just sitting in it.

- Before use, thoroughly read this Operating Manual for proper and safe use.
- Before use, make sure to read "Safety Precautions" .
- Keep this Operating Manual close at hand for easy access at any time.
- This product is planned, developed, designed, assembled, and quality assured at our factory in Japan.
- This product is designed for home use.

Contents

Before Use

| | |
|--------------------------------|----|
| ■ Components & Description | 3 |
| ■ Before Starting Massage | 6 |
| ■ How to Use the Monitor | 11 |
| ■ First-time Personal Settings | 12 |

How to Use the iCOCOON

| | |
|---|----|
| ■ iCOCOON My Position Adjustment Screen | 16 |
| ■ iCOCOON Program Menu | 19 |
| ■ Use the Program. | 27 |
| ■ Sub Button Menu | 30 |
| ■ Bluetooth® Settings | 33 |

| | |
|--------------------|----|
| ■ Care and Storage | 37 |
|--------------------|----|

| | |
|-------|----|
| ■ FAQ | 38 |
|-------|----|

| | |
|-------------------|----|
| ■ Troubleshooting | 41 |
|-------------------|----|

| | |
|-----------|----|
| ■ Service | 43 |
|-----------|----|

| | |
|------------------|----|
| ■ Specifications | 44 |
|------------------|----|

Main Features

Creating the ultimate relaxation space Relaxation through light, music, and 1/f fluctuations

~ An Approach to Sleep Induction and
Awakening through Light and Music ~

Under the supervision of Dr. Obayashi of Nara
Medical University, we have introduced a
light program that is aligned with biological
rhythms which promote sleep and activation.

~ 1/f Fluctuation Feature ~

Our unique technology enables
1/f fluctuation in the seat cushion.
The combination of irregular, delicate
soothing effects and science creates blissful
moments.

~ Gentle Stimulation in the
"Meditation Program" ~

Utilizing massage functions that provide
gentle stimulation, this program relaxes
muscles stiffened by daily fatigue, harmonizes
breathing rhythms, and leads to relaxation of
both body and mind.


Achieves the ultimate in seating comfort
(cushioning).

Posture Support

Eight air cells automatically create a position
that suits you

Components & Description

Main Unit (Front side)



| | |
|---|--------------|
| 1 | Top Panel |
| 2 | Speakers |
| 3 | Pillow |
| 4 | Back Cushion |

| | |
|---|------------------------------|
| 5 | Monitor (with Monitor Cable) |
| 6 | Main Unit Remote Control |
| 7 | Seat Cushion |
| 8 | Seat Heater |

| | |
|----|------------------------------|
| 9 | Foot Heaters(right and left) |
| 10 | Side Panels(right and left) |
| 11 | Foot Slide |
| 12 | Casters (right and left) |

Armrest

| | |
|-----|------------------------------------|
| A-1 | Shoulder Mechanism(right and left) |
| A-2 | Arm Mechanism (right and left) |


Footrest

| | |
|-----|--------------------------------------|
| F-1 | Calf Mechanism (right and left) |
| F-2 | Foot Sole Mechanism (right and left) |

Continued (Components & Description)

Before Use

Main Unit (Back side)



| | |
|----|-----------------------------|
| 13 | Side Panel (right and left) |
| 14 | Caster (right and left) |
| 15 | Power Cord |
| 16 | Main Power Switch |


| | |
|----|------------|
| 17 | Inlet |
| 18 | Connector |
| 19 | Power Plug |

Air Cells


| | |
|-----|---------------------------|
| K-1 | Pillow Air Cell |
| K-2 | Lower Back Right Air Cell |
| K-3 | Lower Back Left Air Cell |
| K-4 | Seat Right Air Cell |






| | |
|-----|-----------------------|
| K-5 | Seat Left Air Cell |
| K-6 | Seat Front Air Cell |
| K-7 | Back Cushion Air Cell |
| K-8 | Seat Cushion Air Cell |

Air Cells



Monitor / Main Unit Remote Control




| | | |
|---|---|---|
| 1 | Monitor (Accessory) | Operate by tapping the screen. |
| 2 | Main Unit Remote Control | You can operate the main unit and also can install the monitor. |
| 3 |  ON/OFF Button | Turns the main unit power "ON/OFF." When the remote control power is "ON," it lights up green. When it is "OFF," the light turns off. When the power turns off, the reclining seat back will rise and massage rollers will return to the storage position before stopping. The reclining seat back automatically rises when pressing the button during the stage of emergency stop. |
| 4 |  Emergency Stop Button | Press this button in case of an emergency or if you sense any abnormality during the operation. The operation will stop immediately. When an emergency stop is being made, the button flashes in red. |
| 5 |  FULL BODY | This is a recommended program that provides a full-body massage that helps with recovery from fatigue. |
| 6 |  Open Button | The side panels will open. They will continue opening while holding down the button. |
| 7 |  Close Button | The side panels will close. They will continue closing while holding down the button. |
| 8 | USB Port | This is a communication port for customer support. Please do not use it. |

Before Starting Massage

Check around the chair

Make sure that no one, pets or objects are around or behind the chair or in front of the leg section.



Turn on the power

Warning

- Use the product only at the specified power source.
→Using the wrong power source may result in a fire or electric shock.
- Do not use the product with a transformer. →This may result in malfunction or electric shock.
- Do not use the product if the power cord or power plug is damaged or its insertion into a wall outlet is loose.
→This may result in an electric shock, short circuit, or fire.

Caution

- When using the product after raising the ambient temperature in a very cold room, do not increase the room temperature suddenly but gradually.
→Condensation may occur on the operating parts, and performance may not be fully demonstrated.


Main Power Connection

- 1 Plug the power cord into the cord connector located at the back of the chair first. Then plug the other end into the wall outlet.



Set the main POWER switch located at the back of the chair to ON.

A text message should appear on the remote control display screen.



* If the main POWER switch is turned on while the rollers are stopped in the middle of the movement, they will automatically return to their storage position (upper part of the reclining seat) then stop. The seat will also move for adjustment of the initial position then stop.

Warning


- Do not use with any type of transformer device.
→May result in damage to the product or electric shock.

Caution


- Before plugging the power cord into the wall outlet, make sure that the main POWER switch is set to OFF.
- Make sure that the plugs on the power cord are inserted fully into the cord connector and the wall outlet.
- Before using the chair for the first time, turn on all the switches in order as indicated in this manual to check that the chair is operating properly.

2 Turn **ON/OFF** the Main Power Switch of the main unit.

- The main power switch is located on the back side of the main unit.



3 The startup screen appears on the monitor screen. During preparation, the LED on the **ON/OFF** button flashes in green.



Check the main unit

1 Make sure to check that the fabric of the main unit is not torn.

Warning

- Before use, make sure to check that the fabric of the main unit is not torn. Also, make sure to check that the fabric of other parts is not torn.
- However small a break, immediately stop using the product, disconnect the power plug, and contact your local distributor. Never remove the fabric of the main unit or use it in a torn state as it is very dangerous.

2 Make sure that no foreign objects are caught in the main unit.

Caution

- Make sure to check that no foreign objects are caught in the main unit before use.
→ May cause an accident or injury.

3 Make sure to check that each mechanism is in its stowed position.

Mechanism stowed positions

- Back mechanism ... Top of the backrest
- Sole mechanism ... Open state

Caution

- Make sure that each mechanism is in the stowed position before sitting down.
→ May cause an accident or injury.

4 Make sure to check that the backrest is fully raised and the footrest is fully lowered.

- * When the power is on, the side panels will open, and air will be supplied to the pillow, back cushion, and lumbar cushion.

Caution

- When the footrest is not fully lowered, do not sit down, stand, or place an object on the footrest.
→ May cause an accident or injury.

Continued (Before Starting Massage)


How to sit down on the chair

Normally, each mechanism is in a stowed state.


Back mechanism ... Bottom of the backrest; Sole mechanism ... Open state

For safety, check each mechanism position before sitting down, and if it is not in a stowed state, do not sit down immediately but turn on the main power switch to return each mechanism to the stowed state, and then slowly sit down in the center of the chair.

The panels are closed



The panels are open





Caution

- Do not sit in the chair while the footrest is raised.
→An impact may be applied to the footrest causing a malfunction. Also, an accident or injury may result in the case of the main unit overturning.
- Before use, fully raise the backrest.

Precautions for getting in/out of the chair

1 Make sure to check that the footrest is fully lowered, and then put your feet into the sole unit.

2 Get in/out of the chair by placing your hands on the upper part of the footrest.



Warning

- Do not subject the side panels, main unit remote control or armrests to strong impacts or lean against them.
→Overturning them may cause an accident, injury, or a malfunction.
- If the arm mechanism and/or footrest is in operation, wait for the operation to stop, and then get in/out of the chair.
- Be sure to sit in the chair with the side panels open.

Caution

- Do not sit down on the footrest while it is raised.
- Never get in/out while each mechanism is moving because it is very dangerous.
→May result in injury or a malfunction.
- We recommend use of this product for persons weighing 100kg (220 lbs.) or less.
When a person weighing 100kg (220 lbs.) sits in the chair, a load of approximately 1.9kg (4.2 lbs.) is applied to the floor per 1cm² (0.16 in²).
- Before use, make sure to check that the fabric is not torn.
- The withstand load of the armrests is approximately 30kg (66 lbs.).
- The tensile load of the armrests is approximately 30kg (66 lbs.).
- The withstand load of the footrest is approximately 30kg (66 lbs.).

Arm positions

Put your arms in between the armrest.

- Do not subject the armrest to strong impacts or lean on it. It may cause deformation or damage.

Correct arm positions



Place your hands palm-down.



Warning

- Do not place your arms outside the elbow rests of the armrests.
- The arm mechanisms slide. Do not put your hands or fingers in the spaces between the arm mechanisms and armrests.
→May result in injury.




Caution

- Do not use the armrests for anything other than your arms.
→May result in injury.
- If one of your arms comes out of the armrest while your arms are being massaged, do not try to force it back in. Undue force may cause a malfunction.
- Do not massage your elbows.
→May result in injury.

Leg positions

1 Put your legs in between the calf mechanisms and footrest as shown in the illustration.

2 Sit back in the seat. It does not matter if your feet are not touching the bottom of the footrest.



Automatic Foot Slide Adjustment Function


The Automatic Foot Slide Adjustment Function extends and retracts the footrest by responding to a sensor in the toe area of the footrest. For sensors to react properly, sit down with your soles firmly pressed against the footrest.

* If your toes do not reach the footrest, make adjustments to an angle of the backrest until your legs fit the footrest. The foot slide will be automatically adjusted when starting the massage program.

Continued (Before Starting Massage)

Extension and retraction of the footrest

When raised to a certain height, the footrest can be extended/retracted.



*The extending/retracting length is approximately 27cm (11 inch).

Warning


- Do not sit on the footrest while it is being raised.
- When lowering the footrest, make sure to check there is no one (child in particular) or any pet around the chair. Also, be careful not to get your body or other objects caught under the footrest.
- Do not suddenly release your hand when lifting the footrest.
- When the foot slide function is in motion, be careful not to get your feet or other objects caught in between the foot slide and the footrest.

Caution

- Do not use the footrest for parts other than your legs.
→May result in injury.
- If your one of your legs comes out of the footrest while your legs are being massaged, do not try to force it back in.
→Undue force may cause a malfunction.
- Do not massage your knees. Do not put your legs in other directions.
→May result in injury.

Monitor Position

Remove the monitor located between the back cushions and attach it to the main unit remote control. The monitor and the main unit remote control are attached together with magnets.




Caution

- Do not use if you have an electronic medical device such as a pacemaker.
- Do not bring any magnetic storage media, such as magnetic cards, near the chair. It may result in destruction of data. Additionally, do not bring PCs, electronic wristwatches, or other precision electronic devices close to the chair, as this may cause them to malfunction.
- Do not move the angle adjustment part of the monitor cable beyond its movable range or try to forcibly pull or bend it.
→ The harness inside the monitor cable may break and cause damage.
- Do not apply any load to the monitor or the main unit remote control, such as leaning against or giving a strong impact to them.
→ This may result in an accident, injury or malfunction.


How to Use the Monitor (Before sitting in the iCOCOON)

How to activate the monitor

1 Turn **ON** the main power switch located at the back left of the iCOCOON main unit.



The startup screen appears on the monitor screen. During preparation, the LED on the **ON/OFF** button flashes in green.




2 When the standby screen appears, press the **ON/OFF** button of the main unit remote control.


- After the elapse of a certain period of time with no operation being performed, the screen turns black. press the **ON/OFF** button to turn it on again.

* For details, refer to the FAQ.

3 The power is turned "ON" and the iCOCOON panels open. When the power is ON, the LED on the **ON/OFF** button lights up green. After that, air is automatically supplied to the pillow, back cushion and seat cushion.



4 The User Selection screen will appear after air preparation is completed. Please sit in the iCOCOON.




First-time Personal Settings


By registering individual users, users can manage their data individually.

User Registration

- 1 Tap the **New Registration** button on the Login screen.
 - You can register up to 8 people.




- 2 Create your account. Tap the **OK** button on the screen. When the **Cancel** button is tapped, the new registration is interrupted and move to the login screen.




- 3 The Name Entry screen will appear. Enter <NickName>. After entry, tap the **Enter** button.




- 4 The Date of Birth Entry screen will appear. Enter year of birth. After entry, tap the **Enter** button.



- 5 Enter month of birth. After entry, tap the **Enter** button.




- 6 Enter day of birth. After entry, tap the **Enter** button.



(Example: For those born in January, please enter '1', and for those born in October, please enter '10'.)


7

Select <Gender>.




8

Select <Height>.




9

Select <Foot size>.




10

Select <Shoulder width>.



11

Check the entered contents. If the contents are correct, tap the **Register** button. When editing the contents, tap the **Edit** button



12


The My Position Adjustment screen will be displayed after completion of user registration.




Continued (First-time Personal Setting)

Deletion of user information

- 1 Tap the **Delete user** button on the Login screen. Select the user you want to delete.




- 2 After confirmation, tap the **Delete** button so that the information is deleted.
 - Note that once erased, the previously registered information cannot be restored.



User Login

- 1 Tap the user button which is registered in the User Selection screen.




- 2 The My Position Adjustment screen will be displayed. Tap the **Program Menu** button in the upper right of the screen to display the Program Menu screen.



- 3 Tap the button in the upper right of the screen to toggle between the Program Menu screen and the My Position Adjustment screen.


iCOCOON Program Menu(Go to Page 20.)



Guest Login

*It is necessary to enter physical information each time. The entered data is not recorded.


- 1 Tap the **Guest** button on the User Selection screen.




- 2 Select **<Foot size>**.



- 3 Select **<Shoulder width>**.




- 4 Select **<Gender>**.



- 5 Select **<Height>**. After selection, tap the **OK** button.




- 6 Once the information has been registered, the My Position Adjustment screen will appear. Guest logins cannot register My Position.




iCOCOON My Position Adjustment Screen

Use as a sofa.


- 1 Tap the registered user button to log in.



- 2 For User Logins, the registered My Position is applied.



- 3 The My Position Adjustment screen will be displayed. You can register the amount of air into the pillow, back cushion and seat cushion as My Position.



To register the iCOCOON seating position as My Position

The air volume for the pillow, back cushion, and seat cushion are set as the default position at level **3**. To adjust each cushion's air volume to your preference, select the desired volume of air for each cushion and tap the **Apply** button. This will adjust the air accordingly for each cushion. Tap the **Register** button if you wish to register your preferred position. Your registered position will be reflected the next time you log in.

*Registration of My Position can only be made for User Logins.

Air Intake button of Air Position



Air volume increases one level with a single button tap. Compared to default level **3**, the sitting comfort gets firm and hard.

Air Release button of Air Position



Air volume decreases one level with a single button tap. Compared to default level **3**, sitting comfort gets soft.

Reclining

The backrest is raised/reclined. Operates only while the button is held down.

Footrest

The footrest is raised/lowered. Operates only while the button is held down.

Foot Slide

Extends/retracts the foot slide. Operates only while the button is held down.

30-Minute Position Change Mode Start button

Tapping the button of the [Position Change Mode Start] every 30 minutes will move to the screen to start the position change mode every 30 minutes, and start the position change mode in which air is automatically supplied to the lower back and the seat. The position change mode lasts for a max. of around 6 hours. After this max. time has elapsed, it will automatically change to the My Position Adjustment screen.


- Movement of the seat air support when the program starts (5 min.)
- ↓
- Stop movement 30 min.
- ↓
- Movement of the lower back air support (6 min.)
- ↓
- Stop movement 30 min.
- ↓
- Movement of the seat air support (5 min.)
- ↓
- Stop movement 30 min.
- ↓
- Then, the series of movements will be repeated.

Release all air

Air continues to be released from the pillow air cell, lower back right & left air cells, back cushion air cell and seat cushion air cell as long as the button is held down.

- *If you feel that there is too much air in the air bags even after adjusting My Position, or if you want to use the chair with even less air than Position 1, use the [Release All Air] button to release the air.
- *If the air is coming out slowly, please press the air cell from above. The applied pressure makes it easier to release the air from the air bags.

- 4 Starts a 30-minute position change mode in which air is automatically supplied to the lower back (lower back left air cell, lower back right air cell) and the seat (front seat air cell, seat left air cell, seat right air cell). During the program, the reclining angle and other features can be adjusted using the buttons on the screen. After the program is completed, the screen automatically changes from the 30-minute Position Change Mode screen to the My Position Adjustment screen including the adjusted reclining angle.



Reclining

The backrest is raised/reclined. Operates only while the button is held down.

Footrest

The footrest is raised/lowered. Operates only while the button is held down.


Foot Slide

The foot slide is extended/retracted. Operates only while the button is held down.


iCOCOON Program Menu

Types of programs


1 Tap the registered user button to log in.



2 Registered My Position is applied.




3 The My Position Adjustment screen will appear. Tap the **Program Menu** button on the upper right of the screen when you use a program.



Continued (iCOCOON Program Menu)

4 The Program Menu screen will appear.



Program Menu is grouped into two types.

①Meditation (3 programs) (Go to Page 21.)

Family Inada's concept of meditation is to relax and harmonize your mind and body. The air gently embraces you, and 1/f fluctuations and other movements are performed.

*"1/f fluctuation" refers to fluctuations in which regularity and irregularity exist in a comfortable balance, and is said to calm the human mind.

②Massage/Stretch (9 programs) (Go to Page 22.)

Massage & Stretching


You can experience the technology that Family Inada has accumulated over the years.

Benefits

As a therapeutic massage device for household use, the iCOCOON provides the following benefits:

- Relaxes muscles
- Relives nerve and muscle pain.

1 Meditation



Three different programs can be selected the purpose for mental and physical health through meditation that makes full use of 1/f fluctuation, light, music, and a massage function that provides gentle stimulation.

Sleep Induction Program

LED

Sleep Induction

After gently massaging the entire body, it embraces you in warm light to induce sleep and lead you to good quality sleep.

Mental Fatigue Recovery Program

LED

Meditation

This program promotes relaxation of muscles that are stiff from daily fatigue, enhances relaxation, and improves the effects of meditation.

1/f Fluctuation Program

LED

Meditation


This program promotes irregular, pleasing 1/f fluctuations to release physical and mental tension, guiding you to deep relaxation.

*"1/f fluctuation" refers to fluctuations in which regularity and irregularity exist in a comfortable balance, and is said to calm the human mind.

Continued (iCOCOON Program Menu)

② Massage/Stretch

How to use the iCOCOON




Full Body Program

LED

Sleep Induction

This program lets you experience the distinctive movements of iCOCOON (exquisite massage and meditation that brings peace of mind and body) together with a unique ambience.



Massage


You can choose from five different programs according to your needs. (Go to Page 23.)

Posture Stretch Program

LED

Wakefulness


This program aims to enhance body flexibility and range of motion, relieve muscle tension, and promote mental and physical relaxation.



Children

You can choose from two different programs for young people aged 10 years and older. (Go to Page 26.)

② Massage/Stretch



Morning Program

LED

Wakefulness

Adopts a technique originated from the ancient European concept of "massage," which enhances blood flow from the extremities of the body towards the heart. It promotes blood circulation, relieves fatigue, increases metabolism, and refreshes the body.

Night Program

LED

Sleep Induction

Adopts a technique originated from the ancient Chinese concept of "Anma," which enhances blood flow from the heart to the extremities of the body. It promotes relaxation.

Stubborn Shoulder Stiffness Program

LED

Wakefulness

This program focuses on the levator scapulae muscle, often blamed for shoulder stiffness, and loosens the tense muscles around the shoulder area. It begins by relaxing the surface muscles of the neck and shoulders, and then applies firm kneading for deeper relief.

Lymph Program


LED

Wakefulness

A massage program that follows the flow of the lymphatic system, utilizing a slow rhythm.

Continued (iCOCOON Program Menu)

② Massage / Stretch



Area Specific Program


This program is a specialized and intensive massage for each area such as legs and the lower back of the body.

*It cannot be combined with other parts of the body.


*IProlonged massage on the same area is not recommended.

*In this program, please select your preferred LED lighting and enjoy choreographed lights.


- 1 Tapping the **Area Specific Program** button automatically reclines the seat back and performs foot sensing and shoulder position detection.



3 Foot Sensing will be performed. As the footrest retracts, please straighten your knees. The foot sole sensor will react automatically, causing the footrest to extend or retract, but you can also adjust it using the buttons on the screen. you select a part, the program for the selected area starts automatically.



4 After foot sensing is completed, the back mechanism rises, and shoulder position detection will be performed. If adjustment of the shoulder position is necessary, use the **Up** **Down** buttons on the screen to adjust the shoulder position.



5 After the shoulder position is detected, a screen for selecting specific body area for Area Specific appears. When you select a part, the program for the selected area starts automatically.

An area you want to unstiffen in a concentrated manner can be massaged. Choose an area you want to massage.



Continued (iCOCOON Program Menu)

② Massage / Stretch



Improve Concentration Program

LED

Wakefulness

This program is incorporates the results of joint research to improve the concentration of young people with Dr. Fukada, Doctor of Medicine, Tottori University School of Medicine.

<Summary of Joint Research>

Comparing before and after massage, the time required to answer math questions was reduced and the number of valid answers increased after massage.

*Based on experimental data from our previous models.

Fatigue Recovery Program

LED

Sleep Induction

Children and adults differ not only in body size but also in skeletal structure and muscle mass. The goal is to provide a soft and comfortable massage to relieve daily fatigue and stress and to provide relief for muscles after exercise, while recognizing the differences in physique.




Caution

- Children's program is designed for children 10 years old or more. Parent must accompany children when using this program. Children below the recommended age must not use iCOCOON.

Use the Program

Massage


1 Select your preferred program button.




How to use the
iCOCOON

2 A program overview screen will be displayed. Tap the preferred program button.


*Depending on the type of program, a program overview screen may not be appeared.



3 A detailed program description screen will be displayed. If you wish to start the program, tap the **Start Program** button.




4 After starting the program, the backrest automatically reclines.




Continued (Use the Program)


- 5 Foot Sensing will be performed.
As the footrest retracts, please straighten your knees. The foot sole sensor will react automatically, causing the footrest to extend or retract, but you can also adjust it using the buttons on the screen.



- 6 After foot sensing is completed, the back mechanism rises, and shoulder position detection will be performed.
If adjustment of the shoulder position is necessary, use the **Up** **Down** buttons on the screen to adjust the shoulder position.



7 Starts the program. During the program, you can adjust the intensity and reclining angle using the buttons on the screen.



Intensity Adjustment

If you feel the massage is too strong or too weak during the program, you can adjust the intensity.

Footrest

The footrest can be raised/lowered. It operates only while the button is held down.

Heater

You can select whether the toe heater/the seat heater is ON / OFF.

Backrest

The backrest can be raised or reclined. It operates only while the button is held down.

Foot Slide

The foot slide can be extended/retracted. It operates only while the button is held down.

Automatic Raise Up

You can select whether or not to automatically raise the reclining backrest after completion of the massage program. If you turn it OFF the backrest does not rise automatically after completion of the massage program.

During a guest login, the Automatic Restore to Upright Position is fixed to ON.

Continued (Use the Program)

Sub Button Menu

Guest

My Position adjustment

Program Menu

Deflate Inflation

Reclining Foot Rest Foot Slide

30-Minute Position Change Mode

Release all air

You can adjust the air volume of iCOCOON.

Side Panel LED (Light Patterns)

The side panel's LED will illuminate according to three different intended light patterns.

Sleep Induction Mode – A mode that promotes sleep.

| | | |
|--|---------------|---|
|  <p>This light pattern is intended to induce sleep.</p> | Light color | Warm Color Light (Orange) |
| | Light Control | Bright to Dark (Gets gradually darker.) |

Wakefulness Mode – A mode that promotes wakefulness.

| | | |
|---|---------------|---|
|  <p>This light pattern is intended to promote wakefulness.</p> | Light color | White Color Light |
| | Light Control | Dark to Bright (Gets gradually brighter.) |

Meditation Mode – A mode that promotes a calm state.

| | | |
|---|---------------|---|
|  <p>This light pattern is intended for meditation.</p> | Light color | Warm Color Light (Red to Orange) |
| | Light Control | Simulates the flickering of a candle flame. |

Sub Button Menu

Side Panels LED Adjustment



The LEDs on the side panels can be illuminated in the color of your choice and the brightness can be adjusted in three levels.

Volume

The volume can be set in 10 levels.



– Turn down the volume

+ Turn up the volume

Back Close the sounddialog.

Mute

Guide

QR CODE for the Operating Manual and the Installation Manual.

Open the Camera app of your device. Hold your Camera so that the QR code appears in the viewfinder in the Camera app. Your device recognize the QR code and shows a notification or link address. Tap it to visit the website of the manuals.




Continued (Use the Program)

Sub Button Menu

Settings


Transition to the various settings screen.

How to use the iCOCOON



Edit User Information

You can edit the registered user information.



Set up Bluetooth®


You can connect your device to the iCOCOON. (For connection setup, refer to Pages 33 to 36 of this Operating Manual.)

Bluetooth® Settings

You can play your favorite music through the iCOCOON speaker by pairing your device with the iCOCOON via Bluetooth®.


* Even with a Bluetooth® connection, you cannot display still images or videos on the monitor.

Connect the iCOCOON and your device (such as a smartphone) via Bluetooth®.

- 1 Connect the iCOCOON and your device via Bluetooth®.
Tap the  button in the sub button menu at the footer.




- 2 Tap the **Bluetooth** button.
The Bluetooth® Settings screen will appear.




- 3 Operate your device to turn on the Bluetooth® feature.
[Example: Smartphone]


Select [Settings].



Select [Bluetooth®].



Turn [Bluetooth®] on.





* Depending on your device, the [Bluetooth®] option may be displayed as [Device Connection] or similar.

* The display content may vary depending on your device.

Continued (Bluetooth® Settings)

4 To register a new device on the iCOCOON, tap the **Device Registration** button. The iCOCOON will switch to Bluetooth® connection mode.


* The Bluetooth® connection mode of the iCOCOON lasts for 60 seconds.




5 * The Bluetooth® connection mode of the iCOCOON lasts for 60 seconds.

[Example: Smartphone]



If a Bluetooth® connection confirmation is displayed, select [Pair Device.]



Select [Pair New Device.] Searching for connectable devices will begin.



A list of connectable devices will be displayed. Please select "iCOCOON_YYYYYY" (Six digits of the serial number).



On your smartphone, once "iCOCOON_YYYYYY" appears as [Registered Device] or [Connected], the registration is complete.


* Depending on your device model, when Bluetooth® is turned on, it automatically searches for connectable devices.

* If no devices are displayed, please search again.

* The six-digit serial number is the same as the number found on the nameplate (affixed to the main unit remote control, see Page 43).

* The display content may vary depending on your device.

6 Once the Bluetooth® connection is complete, the screen will automatically switch, and the name of the registered device will be displayed as [Paired Device] on the Bluetooth® Settings screen.



- Each user can register one device per account. Only one device can be registered per user.
- The registered device will automatically connect to the iCOCOON upon the next login.


Play music stored on your device from the iCOCOON speaker

After confirming that the device name is displayed on the Bluetooth® Settings screen, please play your favorite music or other audio from the device.

[Example: Using a smartphone]

Select your preferred music player app or video app.

Choose your favorite music or video and play it. If the Bluetooth® connection has been established, the music or audio will play through the speaker of the iCOCOON.



- If you want to adjust the volume of the audio streaming from the connected Bluetooth® device, please operate your device.


Continued (Bluetooth® Settings)

Disconnect the Bluetooth® connection between the iCOCOON and the device (such as a smartphone).

Disconnect Button

Tap the **Disconnect** button on the Bluetooth® Settings screen. You can disconnect only the Bluetooth® connection.

The device is registered as [Paired Device] on the Bluetooth® Settings screen.




When reconnecting the iCOCOON and the device, tap the **Connect** button again.

Connection via Bluetooth® is made again automatically.

Delete the registered device (such as a smartphone) from the iCOCOON.

Delete Button

Tap the **Delete** button on the Bluetooth® Settings screen. You can delete the registered device from the Bluetooth® Settings screen.



When reconnecting the iCOCOON and your device, register the device using the

Device Registration button and then connect them via Bluetooth®.

Care and Storage

How to care

Caution

- For care and maintenance, make sure to turn off the main power switch, and unplug the power plug from the wall outlet. Never plug in or unplug the power plug with wet hands.

Genuine leather parts

- Normally, wipe off dust or dirt with a dry cloth.
- Wipe off dirty adherents with a soft cloth wrung out of lukewarm water in which a small amount of neutral detergent has been dissolved.
For the care of genuine leather products, it is recommended to use a cleaner and conditioner specifically designed for leather if possible.
- Tightly wring out a cloth soaked in water and wipe off any remaining detergent from the surface. Use a dry, soft cloth to remove the moisture from the surface, and let it dry in a well-ventilated, shady area.
- To maintain quality for a long period of time, we recommend performing regular maintenance about twice a year.

~Caution~

Please follow these instructions to avoid deterioration and damage to the leather surface.

- Remove any dust or dirt on the leather immediately.
- Refrain from using the product in an area exposed to direct sunlight.
- Do not place vinyl, plastic, or wax-containing items on leather surfaces, as they may adhere to the leather when exposed to high temperatures.

Other plastic parts

- Wipe with a soft cloth wrung out of lukewarm water in which a small amount of neutral detergent has been dissolved, and then wipe thoroughly with a dry cloth so that no detergent remains.
- Do not use alcohol or other solvents as well as polishing powder for care and maintenance because it may result in flaws, discoloration, or cracking.

Monitor parts

Wipe off grime and smudge marks on the tablet with a dry cloth.

- Never wipe with a wet towel, etc. May result in a malfunction.

Storage


Caution

- Do not install the product in a place subjected to high temperature such as an area exposed to direct sunlight or near a stove.


Store the main unit in a place with little dampness and dust after removing dirt and dust. If the product is not used for a long time, cover it to prevent adhesion of dust.

- In order to prevent your pet(s) from biting the cord, do not keep your pet away from this unit during storage.

What is appropriate



What is not appropriate



Q Is it OK if I install the chair on flooring or floor heating?

A There is no problem.

However, since the floor may be damaged, it is recommended to lay a carpet and install the massage chair on it. (The same applies to tatami or the like.)

Do not use an electric carpet because the heating element may be damaged, and ignition may occur.

Q Water was spilled on the product by accident. What should I do?

A An electric shock may occur. Immediately stop using the product, and follow the procedures given below.

- Turn off the main power switch directly behind the main unit.
- Unplug the power plug from the wall outlet. After carrying out the above, consult your local distributor.

Q How can I make the massage stronger?

A Try the following methods.

- Make adjustments with the intensity adjustment button displayed on the screen during the massage.
- Sit back so as to apply your body weight to the backrest.
- Lower the backrest.
Make a stronger massage by the above methods.

Q How can I make the massage weaker?

A Try the following methods.

- Make adjustments with the intensity adjustment button displayed on the screen during the massage.
- Raise the backrest.
- Place a pillow.
Make a weaker massage by the above methods.

Q Is it okay to use this product if I am a heavy person?

A We recommend use of this product for persons weighing 100kg (220lbs) or less to use the chair.

Q Can a person under 150cm (4.92 ft) or taller than 180cm (5.91 ft) in height use the chair?

A Yes.

Since the shoulder position cannot be adjusted at the time of shiatsu point search, We recommend using the following method.

- Short person ... Raise the chair.
- Tall person ... Lower the chair.

Q Is there an age limit for users?

A Small children and people who cannot operate the chair on their own cannot use it for security reasons.

However, only the younger user program (However, only the Children Program can be used by younger user 10 or older. In such case, use by younger users should always be in the presence of a guardian. Elderly persons age 80 or older, who have a high risk of osteoporosis, should consult a doctor before use.

Q How often can I use the massage chair each day?

A As a guide, one person should have a massage (30 minutes once or 15 minutes twice, etc.) each day.

When using Chair Menu other than massage, use exceeding the above is not a problem.

Q The foot slide adjustment function does not work properly. (The footrest does not fit my soles.)

A For sensors at the heel section to react, sit with your heels pressed against the sole unit. If your heels do not reach it, adjustments can be made with the foot slide buttons.

Q The side panels do not move even when pressing the [Open] or [Close] button on the main unit remote control.

A Please check the following.

- The side panel operation is restricted when the iCOCOON performs any of the following actions.
- Power-off state
- Preparing state
- Storage movement in progress
- My Position function in progress

Q Pressing the [Full Body] button on the main unit remote control does not start the massage.

A Please check the following.

- Operation is restricted when the iCOCOON performs any of the following actions
- Power-off state
- Preparing state
- Storage movement in progress
- My Position function in progress

Q What should I do if I want to raise the backrest of the iCOCOON during an emergency stop?

A Press the ON/OFF button to make the backrest automatically rise during an emergency stop.

Q The monitor screen went black automatically even though the main power is on.

A After the elapse of a certain period of time with no operation being performed, the monitor screen automatically turn off certain amount of time. Please check the following

<Power-off State>

- On the standby screen
If there is no operation for about 3 minutes, the monitor backlight will go off. To turn it on again, press the ON/OFF button.

<Power-on State>

- On the My Position Adjustment screen
If there is no operation for about 5 minutes, the monitor backlight will go off. To turn it on again, touch the monitor screen.
- On the 30-minute Position Execution screen
If there is no operation for about 5 minutes, the monitor backlight will go off. To turn it on again, touch the monitor screen.

* Please note that pressing the ON/OFF button in the power-on state will transition to the standby screen, and the main unit's power will turn off.

Q The side panel closed suddenly. What should I do?

A If there is no operation for about one minute on the standby screen after turning the main unit off, the side panel will close. After transitioning to power-off state, the side panel will not close immediately, so take your time getting out of the iCOCOON.

- We recommend using the My Position Adjustment screen when using the iCOCOON as a chair.

Continued (FAQ)

Q When using Bluetooth®, I hear noise on my TV or mobile phone. What should I do?

A Try the following

- Move the noisy device away from the iCOCOON.
- Use different power sources for the iCOCOON and the noisy device.

Q Is it possible to connect multiple Bluetooth® devices?

A You can register one device per user. If a device is already registered, you cannot add a new device. Please delete the registered device and then register a new one.

Q Bluetooth® does not work. What should I do?

A Try the following

- If the distance between the iCOCOON and the Bluetooth® device is too far, pairing may not be successful. Please bring the Bluetooth® device closer to the iCOCOON.
- Ensure that the power is turned on.
 - ① If it is battery-operated, check that the batteries are not out of charge.
 - ② Then, restart the Bluetooth® device.
- Ensure Bluetooth® is ON.
- Please unpair the Bluetooth® device you want to connect and then set up the connection again.
- Is another device connected via Bluetooth® instead of the device you want to connect? If another device is connected, disconnect it and then connect the device you want to connect.

Q Bluetooth® connection was lost. What should I do?

A Devices connected to the iCOCOON via Bluetooth® will automatically disconnect in the following situations. To reconnect, tap the "Connect" button on the Bluetooth® Settings screen.

- The device moves out of Bluetooth® range.
- The power is turned off while the device is connected to the iCOCOON via Bluetooth®.

Q Music via Bluetooth® is not playing from the iCOCOON speaker. What should I do?

A If the device you want to connect to the iCOCOON is already connected to another Bluetooth® device (such as headphones and earphones), the music might not play through the iCOCOON speaker. Disconnect the device from the other Bluetooth® device and then connect it to the iCOCOON via Bluetooth® again.


Q To what degree is extra cost required at the time of a malfunction?

A Make an inquiry to your local distributor.

In the case of abnormal use or excessive use, extra cost may be incurred.

Troubleshooting

Check the following items, and if there still are abnormalities, immediately stop use of the product to prevent accidents, unplug the power plug from the wall outlet, and consult with your local distributor.



Continued (Troubleshooting)

Sound during operation

The following noise types will be heard during use, but these are unavoidable due to the structure and do not affect performance.

- Sound during tapping motion ▶ (Rattling/clattering sound)
- Sound of rollers rubbing against the seat cover during kneading motion ▶ (Squeaking/grinding sound)
- Sound of rollers moving up and down/Sound during kneading motion ▶ (Whirring/rolling/knocking sound)
- Sound during air massage ▶ (Tapping sound)
- Sound of air pump ▶ (Booming sound)
- Sound of air discharge ▶ (Hissing sound)
- Sound of the backrest moving up and down ▶ (Whirring sound)

* These types of sound will not affect the performance, etc.



Warning

- In the case of a malfunction of the equipment, do not repair by yourself, but contact your local distributor.
→May result in a fire, electric shock, or injury.



Caution

- If no movements are made or you sense something abnormal, immediately stop use of the product to prevent an accident, turn off the main power switch, unplug the power plug from the wall outlet, and ask your local distributor for inspection and repair.
→May result in an electric shock or fire due to current leakage, short circuit, etc.

Warranty and After-sales Service

When asking for repair


Before asking for repair, carefully read the section “Troubleshooting” on P41. After checking, if there still are abnormalities, contact/ask your local distributor.

Note that a charge may be incurred for maintenance service.

- (1) Name, address, and phone number
- (2) Serial number
See the company nameplate affixed the Main Unit Remote Control.
- (3) Date of purchase (or date of installation)
- (4) Malfunction condition (Please specify as detailed as possible.)
(Example) Please inform us of the code number indicated on the tablet.

Location of company nameplate

The six-digit serial number
company nameplate



During the warranty period

If you have any questions about after-sales service, contact your local distributor.

After expiration of the warranty period

Consult your local distributor. If product functions can be maintained by repair, repair will be provided at a charge upon your request.

Specifications

| | | | | |
|---------------------------------|----------------|--|--|---------------|
| Product name | | HCP-ICO10000A | HCP-ICO10000D | HCP-ICO10000E |
| Power Supply | | AC 120V | AC 220-240V | AC 110V |
| Rated Power Consumption | | 270W (with the heater used); heater only: 30W | | |
| Rated Frequency | | 60 Hz | 50-60 Hz | 60 Hz |
| Rated Time | | 30 minutes | | |
| Appliance Class | | Class 0I equipment | | |
| Upper Body | Back Mechanism | Kneading Speed | Approx. 5 times/min. at minimum, approx. 20 times/min. at maximum | |
| | | Tapping Speed | Approx. 165 times/min. at minimum, approx. 450 times/min. at maximum | |
| Automatic Program | | 11 types | | |
| Specific Program | | 7 types | | |
| Air Pressure | | Approx. 35kPa | | |
| Footrest Adjustment Range | | Approx. 27cm (11 inch) | | |
| Reclining Angle | | Approx. 110° to approx. 164° from floor level, stepless (electric reclining) | | |
| Dimensions of Main Unit | | Approx. 95(125) × 155 (240) × 150 (95) cm (W/D/H) (Approx. 37(49) × 61 (94) × 59 (37) inch (W/D/H)) Regarding depth and height, numbers in brackets apply when the chair is fully reclined and the footrest is horizontal. | | |
| Weight | | Approx. 145kg (319 lbs) | | |
| Outer covering | | Genuine leather | | |
| Accessories | | Power cord, T-wrench, Caster stopper×2pcs | | |
| Packaging | | 2 packs/chair | | |
| Manufactured and Distributed by | | FAMILY INADA CO., LTD. 1220, Takata, Daisen-cho, Saihaku-gun, Tottori | | |
| Manufactured by | | FAMILY INADA CO., LTD. | | |

MEMO

MEMO



<https://www.family-chair.co.jp>

Head office: 4-2-10, miyahara, Yodogawa-ku, Osaka City, Osaka

Manufactured and Distributed by: FAMILY INADA CO., LTD. 1220, Takata, Daisen-cho, Saihaku-gun, Tottori
(Daisen Factory)

* "Specifications are subject to change without notice due to product improvements
(including designs of the tablet screens)."

DATE OF ISSUE (day/mo/yr) : 06/11/2024(first edition)